

Article

A TaqMan[®] Assay Allows an Accurate Detection and Quantification of *Fusarium* spp., the Causal Agents of Tomato Wilt and Rot Diseases

Maria Doroteia Campos ^{1,*} , Carla Varanda ¹ , Mariana Patanita ¹ , Joana Amaro Ribeiro ¹ , Catarina Campos ¹ , Patrick Materatski ¹ , André Albuquerque ¹  and Maria do Rosário Félix ²

¹ MED—Mediterranean Institute for Agriculture, Environment and Development & CHANGE—Global Change and Sustainability Institute, Institute for Advanced Studies and Research, Universidade de Évora, Pólo da Mitra, Ap. 94, 7006-554 Évora, Portugal

² MED—Mediterranean Institute for Agriculture, Environment and Development & CHANGE—Global Change and Sustainability Institute, Departamento de Fitotecnia, Escola de Ciências e Tecnologia, Universidade de Évora, Pólo da Mitra, Ap. 94, 7006-554 Évora, Portugal

* Correspondence: mdcc@uevora.pt

Simple Summary: In tomato plants, *Fusarium* spp. have been increasingly associated with several wilt and rot diseases that are responsible for severe yield losses. In this context, a molecular-based tool was developed to increase the accuracy of detection and quantification of *Fusarium* spp. genomic DNA (gDNA) in tomato plants and to discriminate *Fusarium* spp. from other fungal species that affect tomato. This assay revealed to be highly specific and sensitive for *Fusarium* species. The used methodology also allowed the establishment of an absolute DNA quantification method. Finally, the effectiveness of the assay was successfully validated with the detection and quantification of *Fusarium* spp. in potentially infected tomato plants from an experimental field and in control plants grown under controlled conditions. The established methodology allows a reliable, sensitive, and reproducible estimation of *Fusarium* accumulation in infected tomato plants, gaining new insights for disease control and providing an additional tool in the screening of resistant plants.



Citation: Campos, M.D.; Varanda, C.; Patanita, M.; Amaro Ribeiro, J.; Campos, C.; Materatski, P.; Albuquerque, A.; Félix, M.d.R. A TaqMan[®] Assay Allows an Accurate Detection and Quantification of *Fusarium* spp., the Causal Agents of Tomato Wilt and Rot Diseases. *Biology* **2023**, *12*, 268. <https://doi.org/10.3390/biology12020268>

Academic Editor: Juan L. Ramos

Received: 9 January 2023

Revised: 4 February 2023

Accepted: 6 February 2023

Published: 8 February 2023



Copyright: © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<https://creativecommons.org/licenses/by/4.0/>).

Abstract: In tomato plants, *Fusarium* spp. have been increasingly associated with several wilt and rot diseases that are responsible for severe yield losses. Here, we present a real-time PCR TaqMan[®] MGB (Minor Groove Binder) assay to detect and discriminate *Fusarium* spp. from other fungal species that affect tomato plants. The methodology used is based on the selective amplification of the internal transcribed spacer (ITS) region of *Fusarium* spp. This assay revealed to be highly specific and sensitive for *Fusarium* species, targeting only the 29 *Fusarium* isolates from the 45 tested isolates associated to tomato diseases. Sensitivity was assessed with serial dilutions of *Fusarium* genomic DNA, with the limit of detection of 3.05 pg. An absolute DNA quantification method was also established, based on the determination of the absolute number of target copies. Finally, the effectiveness of the assay was successfully validated with the detection and quantification of *Fusarium* spp. in potentially infected tomato plants from an experimental field and in control plants grown under controlled conditions. The established methodology allows a reliable, sensitive, and reproducible estimation of *Fusarium* accumulation in infected tomato plants, gaining new insights for disease control and providing an additional tool in the screening of resistant plants.

Keywords: *Solanum lycopersicum*; *Fusarium* diseases detection; TaqMan[®] probe assay; qPCR

1. Introduction

Tomatoes (*Solanum lycopersicum*), with 187 Mt of production and 5.05 Mha of cultivated area (<http://www.fao.org/faostat/en/>, accessed on 25 November 2022), are the second most cultivated vegetable crop throughout the world following potatoes [1]. However,

they are affected by different pathogens such as fungi, viruses, viroids, oomycetes, bacteria, and nematodes that reduce yield and affect product quality. These pathogen agents cause wilting, leaf spots/blight, fruit spots, and rotting [2].

Among the fungal diseases that affect tomato, we highlight the ones caused by *Fusarium* species. *Fusarium* spp. are ubiquitous soil-borne fungal pathogens which cause destructive vascular wilting, rotting, and damping-off diseases [3–5]. As soil-borne pathogens, *Fusarium* species produce persistent resting structures that are able to survive in the soil for long periods of time in the absence of host plants as dormant propagules (chlamydospores), with their germination triggered by the presence of the host's roots [3,6]. This feature makes eradication quite difficult since chlamydospores may serve as a source of inoculum for infection in subsequent growing seasons [7]. Once the mycelium invades the root cortical cells, the fungi tend to colonize exclusively inside the vessels of xylem, rapidly colonizing the host and producing microconidia, which are transported upwards through the sap stream upon detachment. The characteristic wilt symptoms appear due to vessel blockage triggered by the gathering of fungal hyphae and a combination of the host and pathogen [3,6]. In addition to the losses caused before or during harvest, some *Fusarium* species have the ability to produce mycotoxins that are increasingly becoming a public health concern due to the frequent contamination of various food commodities [8,9].

Fusarium oxysporum is an important fungus, known for its phylogenetic diversity [10]. Considering the *F. oxysporum* species, we highlight the two most important formae speciales *F. oxysporum* f. sp. *lycopersici* (FOL) responsible for Fusarium wilt and *F. oxysporum* f. sp. *radicis-lycopersici* (FORL) responsible for Fusarium crown and root rot diseases. Both formae speciales cause extensive production losses in tomato fields and greenhouses and continue to present major challenges for tomato production [3,11]. In addition to *F. oxysporum*, other *Fusarium* species are also known to cause extensive worldwide crop losses such as *Fusarium equiseti*, *Fusarium graminearum*, *Fusarium proliferatum*, *Fusarium solani*, and *Fusarium verticillioides*, all of them associated with tomato wilt [12–16].

The need to increase food production together with the demand to reduce the application of synthetic chemicals that have consequences on increasing the carbon footprint and negative impacts on the environment and human health, led to the search for alternative methods to protect plants against pathogens. Appropriate management practices are invaluable in reducing plant disease losses and have already demonstrated effectiveness in reducing the presence of *Fusarium* species [11,17]. Biological control using antagonistic microorganisms may also constitute an alternative treatment to control plant diseases with quite promising prospects in several plant species, including tomato [10,18–21]. Moreover, the identification of candidate genes in susceptible and resistant responses may facilitate genetic engineering efforts to incorporate new sources of resistance in tomato plants for protection against pathogens for sustainable plant disease management [4,16]. The use of the highly sensitive real-time quantitative PCR (qPCR) arises as an extremely useful tool for studying various agents of infection in plants, such as fungi, viruses, or bacteria [22–24], leading to a better control of diseases and limiting the use of chemical defence strategies.

Given the high incidence of diseases in tomato plants caused by *Fusarium* spp., their consequent negative economic impacts, and the fact that most phytosanitary treatments are based on the application of synthetic fungicides, the establishment of a molecular-based tool that enables their early and accurate detection is of great interest. Furthermore, it will provide an additional tool for the screening of resistant plants. Until now, qPCR has been applied for the detection and quantification of *Fusarium* species individually or targeting few species at the same time [7,25,26]. In the study presented here, a new TaqMan®-based qPCR method targeting the *Fusarium* spp.-specific internal transcribed spacer (ITS) region was developed for the simultaneous detection and quantification of a panoply of *Fusarium* species that affect tomato plants. As a proof of principle, the new qPCR assay was used to assess *Fusarium* spp. contamination of tomato field plants and of plants grown under controlled conditions. qPCR combined with the chemistry of TaqMan®

MGB probes represents a highly specific and sensitive detection system, even when low amounts of target DNA are present, as in the case of early plant–fungi interactions.

2. Materials and Methods

2.1. Fungal Isolates

For specificity testing in qPCR assays, potato dextrose agar (PDA) plugs of 29 *Fusarium* spp. isolates and 17 non-*Fusarium* isolates were retrieved from long-term storage at 4 °C, transferred to PDA plates, and maintained at room temperature prior to gDNA extraction. The fungal isolates used in this study belong to the collection of the Mycology Laboratory, Mediterranean Institute for Agriculture, Environment and Development (MED), University of Évora, Portugal, and are listed in Table 1.

Table 1. List of isolates of *Fusarium* spp. and other fungal isolates used to test species specificity using the *Fusarium* spp.-specific qPCR assay. Data are expressed as quantification cycle (Cq) values, using 100 ng of gDNA as template. N.D.: not detected.

Species	Isolate Ref.	Cq Value
<i>F. incarnatum</i>	M_4F	22.90
<i>F. oxysporum</i>	M_45-2	22.86
<i>F. oxysporum</i>	M_67	20.48
<i>F. oxysporum</i>	M_84	21.53
<i>F. oxysporum</i>	M_90	25.05
<i>F. oxysporum</i>	M_91	20.31
<i>F. oxysporum</i>	M_92	20.26
<i>F. oxysporum</i>	M_100	19.29
<i>F. oxysporum</i>	A_56	25.08
<i>F. oxysporum</i>	A_57	22.76
<i>F. oxysporum</i>	A_122	23.57
<i>F. oxysporum</i>	A_127a	26.42
<i>F. oxysporum</i>	A_127b	30.32
<i>F. oxysporum</i>	A_137	25.65
<i>F. oxysporum</i>	C_4	18.42
<i>F. oxysporum</i> f.sp. <i>radicis-lycopersici</i>	C_3	19.20
<i>F. solani</i>	M_215c	21.51
<i>F. verticillioides</i>	A_118	29.29
<i>Fusarium nelsonii</i>	A_82c	20.17
<i>Fusarium nelsonii</i>	A_82e	20.17
<i>Fusarium</i> sp.	M_62	20.01
<i>Fusarium</i> sp.	M_207	20.16
<i>Fusarium</i> sp.	M_208	22.19
<i>Fusarium</i> sp.	M_209	20.28
<i>Fusarium</i> sp.	P_Q7	19.95
<i>Fusarium</i> sp.	A_123	24.63
<i>Fusarium</i> sp.	A_126	25.91
<i>Fusarium</i> sp.	C_1	21.48
<i>Fusarium</i> sp.	C_2	21.06
<i>Alternaria alternata</i>	A_59	N.D.
<i>Alternaria alternata</i>	A_72	N.D.
<i>Alternaria alternata</i>	C_6	N.D.
<i>Alternaria tenuissima</i>	P_75B	N.D.
<i>Botrytis cinerea</i>	A_43	N.D.
<i>Botrytis cinerea</i>	A_64	N.D.
<i>Botrytis cinerea</i>	A_115	N.D.
<i>Botrytis cinerea</i>	A_119	N.D.
<i>Botrytis cinerea</i>	A_132	N.D.
<i>Botrytis cinerea</i>	C_9	N.D.
<i>Cladosporium cladosporioides</i>	P_70D	N.D.
<i>Colletotrichum</i> sp.	A_71	N.D.
<i>Epicoccum nigrum</i>	P_R46	N.D.
<i>Phytium</i> sp.	C_10	N.D.
<i>Phytophthora</i> sp.	C_5	N.D.
<i>Verticillium dahliae</i>	C_8	N.D.

2.2. Plant Sampling

The tomato plants used in all experiments belong to the UG 29814 variety from United Genetics (Parma, Italy) and have high and intermediate resistance to several diseases, including *Fusarium* wilt caused by FOL race 3.

For *Fusarium* spp. quantification in tomato tissues, plants from a commercial nursery were planted in a field located in the Ribatejo region (central Portugal) (39.0519608° N, 8.7843653° W). This field has been intensively cropped with tomato plants and known to be infested with *Fusarium* diseases for many years. At the same time, control tomato plants were planted in plastic pots in a mixture of sand and vermiculite (3:1 ratio), autoclaved, and maintained in a room under controlled temperature (22–25 °C) with a 14 h photoperiod.

A total of 30 tomato plants were randomly collected from the experimental field in mid-July, 12 weeks after plantation. At the same timepoint, 16 tomato control plants were also collected. Plant crowns from all samples were detached and the surfaces were disinfected as previously described [27], ground into powder with liquid nitrogen, and stored at −80 °C until further processing.

2.3. gDNA Extraction from Fungal Isolates and Tomato Samples

For each fungal isolate, approximately 100 mg of mycelium from a 2- to 3-week-old PDA culture plate was used for gDNA extraction. gDNA was extracted using the CTAB (hexadecyltrimethylammonium bromide) method [28,29]. For tomato samples, gDNA extraction was performed using DNeasy® Plant Pro Kit (Qiagen, Hilden, Germany) following the manufacturer's instructions. PCR inhibitors from tomato plants' gDNA were removed by the OneStep™ PCR Inhibitor Removal Kit (Zymo Research, Freiburg im Breisgau, Germany), also according to the manufacturer's instructions.

To quantify and assess gDNA purity, the absorbance was evaluated in a Quawell Q9000 micro spectrophotometer (Quawell Technology, Beijing, China). All samples were diluted to a final concentration of 20 ng µL^{−1}. gDNA integrity was checked by 0.8% agarose gel electrophoresis.

2.4. Design of a *Fusarium* spp.-Specific qPCR Assay

Partial sequences of the ITS region from *Fusarium* spp. (including *F. acuminatum*, *F. clavum*, *F. delfinoides*, *F. equiseti*, *F. graminearum*, *F. incarnatum*, *F. oxysporum* f. sp. *cubense*, *F. oxysporum* f. sp. *lycopersici*, *F. oxysporum* f. sp. *radicis-lycopersici*, *F. proliferatum*, *F. sacchari*, *F. solani*, *F. subglutinans*, and *F. verticillioides*) (Figure 1), mostly isolated from tomato plants, were retrieved from the National Centre for Biotechnology Information (NCBI) (<http://www.ncbi.nlm.nih.gov/> accessed on 6 January 2023). Sequences were aligned using the MUSCLE method [30] integrated in the CLC Genomics Workbench 11 (Qiagen, Hilden, Germany). Nucleotide sequences that showed a specific consensus for *Fusarium* spp. were chosen and used with Primer Express 3.0 software (Applied Biosystems, Foster City, CA, USA), selecting the option MGB TaqMan® probes and the default parameters of the software to design a specific probe (FusProbe: 5'-GTTGCCTCGGCGG-3') (Figure 1). The 5' end of the *Fusarium* TaqMan® probe was labelled with 6-FAM (fluorescein), bearing at the 3' end a nonfluorescent quencher (NFQ) coupled to the Minor Groove Binder moiety (MGB) (Eurogentec, Seraing, Belgium). ITS1 and ITS2 primers (ITS1: 5'-TCCGTAGGTGAACCTGCGG-3'; ITS2: 5'-GCTGCGTTCTTCATCGATGC-3') [31] were selected as the probe's flanking primers, generating an amplicon size of 227 bp. To ensure the specificity of the *Fusarium* spp. assay and to demonstrate the suitability of the ITS1-ITS2 region for the design of the specific assay, a bioinformatic analysis was performed, using the correspondent target partial sequence from the ITS region of the *Fusarium* species, as well as from other fungal species previously reported to affect tomato plants in the Mediterranean Basin [1]. A Maximum Likelihood phylogenetic analysis based on the Tamura–Nei model of MEGA 11 software [32] was performed using a bootstrap of 1000 replicates.




Figure 1. Multiple alignment of a partial sequence of the ribosomal internal transcribed spacer (ITS) region of nuclear rDNA of *Fusarium* species and respective GenBank accession numbers. The location of the probe is indicated with the arrow and corresponds to the position 80 bp from the complete sequence of *Fusarium oxysporum* f. sp. *radicis-lycopersici* (accession number MH865886.1).

2.5. Specificity, Sensitivity, and Reliability of the qPCR Assay

All qPCR reactions were performed in 96-well plates in a Line Gene 9600 Plus real-time PCR detection system (BIOER, Hangzhou, China). Reaction mixes contained 100 ng of gDNA as the template, 2× SensiFAST Probe Hi-ROX Kit (Meridian Bioscience, Newtown, OH, USA), 400 nM of each primer, and 100 nM of the probe (Eurogentec, Liège, Belgium) in a total volume of 20 µL. The quantification cycle (Cq) values were obtained for each sample with the following cycling conditions: 20 s at 95 °C for initial denaturation, an amplification program of 40 cycles at 95 °C for 15 s, and 60 °C for 20 s. The fluorescence threshold was manually set to 150. Three technical replicates were considered for each sample. *Fusarium* sp. positive controls and no template controls were included in all plates.

For species specificity validation, a *Fusarium* assay was performed on the gDNA of the 45 described isolates, most of them commonly associated with tomato diseases, including 29 *Fusarium* spp. isolates (Table 1). As a measure of sensitivity and the quantitative range of the developed qPCR procedure, the limit of detection (LOD) was determined. A total of 20 standards were prepared by a two-fold serial dilution of the gDNA of *F. oxysporum* f. sp. *radicis-lycopersici* (2^{-1} , 2^{-2} , ..., 2^{-20} , starting with 100 ng) in a background of tomato gDNA from a healthy plant (Table 2). To determine the reliability of the assay, the detection and quantification of *Fusarium* spp. was performed in 30 potentially infested tomato plants from the experimental field and in 16 control plants (for plant sampling details, see Section 2.2) (Table S1).

2.6. gDNA Calibrator Plasmid

The method used for the absolute DNA quantification was based on the determination of the absolute number of target copies (TCN) present in each sample, corresponding to haploid genome equivalents [33]. To obtain a specific calibrator plasmid to be used on qPCR, a gDNA region of *F. oxysporum* f. sp. *radicis-lycopersici* that comprised the corresponding TaqMan® target sequence was amplified using the ITS1 and ITS4 primers (ITS1 sequence, see 2.4; ITS4: 5'-TCCTCCGCTTATTGATATGC-3') [31] by end-point PCR, using the DreamTaq™ DNA polymerase (Thermo Scientific, Waltham, MA, USA) as previously described [34]. The separation of the generated amplicons was performed by agarose gel electrophoresis, and the amplicon with the size of 545 bp was cloned into a pGem®-T Easy vector (Promega, Madison, WI, USA) and used to transform *Escherichia coli* JM109 (Promega, Madison, WI, USA) competent cells by standard procedures. Plasmid DNA extraction was performed using GeneJET Plasmid Miniprep Kit (Thermo Scientific, Waltham, MA, USA), and the screening of the selected putative recombinant clones was conducted using EcoRI restriction enzyme (Thermo Scientific, Waltham, MA, USA). Plasmid DNA quantification was determined in a Quawell Q9000 micro spectrophotometer (Quawell Technology, Beijing, China). The sequence of the insert of selected recombinant clones was

confirmed by Sanger sequencing (Macrogen, Inc., Madrid, Spain: www.macrogen.com) using the T7 and SP6 universal primers located at the cloning vector (flanking the insert at 5' and 3'). One bacterial clone harbouring each of the recombinant plasmids was chosen.

Table 2. Sensitivity of the *Fusarium* spp.-specific qPCR assay using a serial dilution of gDNA from *F. oxysporum* f. sp. *radicis-lycopersici*. The lowest amount of detection for both assays can be observed. Data are expressed as quantification cycle (Cq) values and target copy number (TCN). P.C.: positive control; SD: standard deviation; N.D.: not detected.

Dilution	gDNA in PCR (ng)	Cq Value (\pm SD)	TCN
P.C.	100.00	14.67 (\pm 0.22)	39,606,399.4
2 ⁻¹	50.00	15.53 (\pm 0.27)	22,779,380.0
2 ⁻²	25.00	16.72 (\pm 0.27)	10,595,944.3
2 ⁻³	12.50	17.52 (\pm 0.29)	6,333,969.6
2 ⁻⁴	6.25	18.57 (\pm 0.27)	3,223,887.6
2 ⁻⁵	3.13	19.61 (\pm 0.17)	1,651,494.4
2 ⁻⁶	1.56	20.50 (\pm 0.08)	931,695.0
2 ⁻⁷	7.81×10^{-1}	21.73 (\pm 0.26)	422,375.1
2 ⁻⁸	3.91×10^{-1}	23.09 (\pm 0.13)	176,120.6
2 ⁻⁹	1.95×10^{-1}	24.52 (\pm 0.25)	70,205.1
2 ⁻¹⁰	9.77×10^{-2}	25.82 (\pm 0.19)	30,425.7
2 ⁻¹¹	4.88×10^{-2}	26.74 (\pm 0.07)	16,836.7
2 ⁻¹²	2.44×10^{-2}	28.38 (\pm 0.16)	5863.5
2 ⁻¹³	1.22×10^{-2}	30.27 (\pm 0.34)	1738.7
2 ⁻¹⁴	6.10×10^{-3}	32.52 (\pm 0.59)	409.0
2 ⁻¹⁵	3.05×10^{-3}	33.42 (\pm 0.29)	229.3
2 ⁻¹⁶	1.53×10^{-3}	N.D.	N.D.

A ten-fold dilution series of the selected recombinant plasmid DNA was used to draw a six-point calibration curve in the dynamic range chosen (8×10^2 to 8×10^7 TCN), using the qPCR conditions as described above. qPCR amplification efficiency was calculated using the formula $E = (10^{(-1/\text{slope})} - 1) \times 100$ [35], as well as slope and linearity (regression coefficients, R^2). *Fusarium* spp. gDNA quantification was carried out by interpolating the samples' Cq values onto the standard curve.

3. Results

3.1. Specificity and Sensitivity of the *Fusarium* spp.-Specific qPCR TaqMan® Assay

A specific *Fusarium* qPCR assay that specifically targets the ITS region of *Fusarium* spp. in tomato plants was successfully developed, allowing the detection and discrimination of *Fusarium* spp. gDNA from other fungal isolates. The specificity was first demonstrated in silico, performing searches against the NCBI database with the probe targeting *F. oxysporum* f. sp. *radicis-lycopersici*, *F. oxysporum* f. sp. *lycopersici*, *F. solani*, *F. oxysporum* f. sp. *cubense*, *F. incarnatum*, *F. equiseti*, *F. graminearum*, *F. verticillioides*, *F. subglutinans*, *F. proliferatum*, *F. sacchari*, *F. clavum*, *F. delfinoide*, and *F. acuminatum* (Figure 1). The phylogenetic analysis that included *Fusarium* spp. and other fungal species affecting tomato plants in the Mediterranean Basin revealed that the ITS1-ITS2 region forms a *Fusarium* spp. cluster and separates it from the other fungal species (Figures 2 and S1).




Figure 2. Phylogenetic tree of a partial sequence of the ribosomal internal transcribed spacer (ITS1–ITS2) region of nuclear rDNA by Maximum Likelihood of *Fusarium* spp. and other fungal species that affect tomato plants with respective GenBank accession numbers.

To ensure its specificity, the assay was then evaluated experimentally using 16 non-*Fusarium* spp. and 29 *Fusarium* spp. isolates. All *Fusarium* isolates were amplified by the assay, with Cq values varying from 18.42 to 30.32, while no amplification was detected in the remaining isolates (Table 1).

When the sensitivity of the assay (expressed by the LOD) was evaluated, the lowest concentration of gDNA detected corresponded to 3.05 pg of *F. oxysporum* f. sp. *radices-lycopersici* with a Cq value of 33.42 (± 0.29) (Table 2), achieved in the 2^{-15} dilution. In the remaining points of the dilution series (from 2^{-16} to 2^{-20}), no amplification was detected. The LOD, which corresponds to the concentration that can be detected with reasonable certainty (95% probability) in the analytical procedure [36], is defined in qPCR as the spike amount of the target organism in dilution that could be detected in more than 95% of replicates [37]. No interference of host gDNA was observed when a healthy tomato plant was used in the dilutions instead of water (result not shown). The dilution series standard curve showed a linear correlation ($R^2 = 0.993$) between Cq and the amount of

gDNA template (Figure S2). This study was performed using *F. oxysporum* gDNA due to the prevalence that this species has among the *Fusarium* species that affect tomato plants.

3.2. Calibration Curves for Quantification of *Fusarium* spp. gDNA

A calibration curve of the 10-fold dilution series of the selected recombinant plasmid DNA was established using the *Fusarium* spp. qPCR assay, which enabled the detection of 8×10^2 TCN, corresponding to a Cq value of 31.73. The calibration curve was characterized by the following parameters: slope (-3.58), Y-intercept (41.87), PCR efficiency (90.2%), and linearity ($R^2 = 0.999$) (Figure 3).




Figure 3. Calibration curve obtained with the *Fusarium* spp.-specific qPCR TaqMan[®] assay showing copy number of template vs. quantification cycle (Cq). (A) Calibration curve constructed based on Cq values obtained from the log of ten-fold dilution series of target plasmid DNA in the dynamic range of 8×10^2 to 8×10^7 target copies. (B) Amplification plot where the cycle number increases as target plasmid DNA concentration decreases. R^2 : regression coefficient.

The equation of the calibration curve was used to infer the gDNA amount of *Fusarium* species in terms of TCN. In the evaluation of the sensitivity of the assay, the gDNA amount in terms of TCN was then determined in each point of the dilution series (Table 2). A TCN of 229.3 corresponding to 3.05 pg of *F. oxysporum* f. sp. *radicis-lycopersici* was the lowest amount detected by the qPCR assay (Table 2).

3.3. Applicability of the *Fusarium* spp.-Specific qPCR TaqMan[®] Assay in Tomato Plants

As a proof of concept to evaluate the practical robustness and accuracy of the assay, tomato plants from the experimental field and grown under controlled conditions were sampled, and the *Fusarium* gDNA amount was evaluated. Our results allowed the detection and quantification of *Fusarium* spp. in twelve of the thirty tested samples from the experimental field and in eight of the sixteen samples of the tomato plants grown under controlled conditions (Table S1).

4. Discussion

The precise detection and identification of plant-infecting fungi is essential to facilitate effective management of diseases, with DNA-based methods allowing an accurate plant disease diagnosis [38]. Furthermore, when combined with information on plant disease severity, the quantification of fungi offers an additional tool in the screening of resistant plants to fungal diseases [26,39].

In the work presented here, a new *Fusarium* spp.-specific qPCR TaqMan[®] assay that targets the ITS region was established for the detection and quantification of *Fusarium* spp. in tomato plants, following procedures previously defined for the use of qPCR in microbial diagnostics [37,40]. qPCR combined with the chemistry of TaqMan[®] MGB probes represents the most specific and sensitive detection system, even in the presence of few amounts of target gDNA [33].

Although the ITS region is quite conserved among fungi, the phylogenetic analysis revealed that the ITS1-ITS2 region was appropriate for the design of the assay since a *Fusarium* spp.-specific cluster was identified (Figure 2), and in silico analysis confirmed the specificity of the assay for *Fusarium* species that might affect tomato plants. In fact, the use of this region for the design of assays benefits from the large number of available sequences in databases, which several studies successfully established [26,41,42]. The specificity of the assay was also verified experimentally using a considerably large group of *Fusarium* and non-*Fusarium* isolates (Table 1). However, a range of Cq values varying from 18.42 to 30.32 was observed among the different *Fusarium* isolates tested (although in 20 of the 29 *Fusarium* isolates, Cq values varied from 18 to 23). This finding might be explained by the existence of single-nucleotide polymorphisms in the primers target region that might affect the efficiency of the reaction. It is also possible that the number of rDNA repeats influences the sensitivity of the detection [43]. As further stated by Lavrinienko et al. [44], the rDNA copy number is not necessarily a species-level trait and it can show great intraspecific variability, but still, this locus is not redundant regarding communities' composition and quantification, requiring a careful interpretation of the results.

The importance of the specificity of the assay is even more relevant in plants growing in natural field conditions that might be infected by multiple pathogens. qPCR has been commonly used for the detection and quantification of *Fusarium* in plants, but often targeting individual *Fusarium* species or few *Fusarium* species at the same time [7,16,25,26,45,46]. However, a single qPCR assay for the broad number of *Fusarium* species that might infect tomato plants, lacking until now, is here applied to this specific group of fungi.

When a calibration curve for the quantification of the *Fusarium* spp. gDNA amount was established, the newly designed assay revealed that E, R², and slope were consistent with the acceptance criteria [47], confirming the accuracy and linear response of the assay over a wide range of dilutions and suggesting the absence of PCR inhibitors. A critical step for an efficient gDNA amplification is to eliminate PCR inhibitors such as the high levels of phenolic compounds present in tomato plants [48] that might affect the sensitivity of the assay [23,49]. PCR inhibitors might be difficult to separate from gDNA during the extraction procedure and can result in false negatives or inaccurate quantification [25]. So, in the studies of the sensitivity of the assay in which the *F. oxysporum* f. sp. *radicis-lycopersici* was diluted in tomato gDNA, as well as in the studies on the applicability of the assay to tomato plants, the use of the DNeasy® Plant Pro Kit combined with the OneStep™ PCR Inhibitor Removal Kit was essential for tomato gDNA extraction.

Although FOL and FORL are among the most intensively studied *Fusarium* species that affect tomato plants, responsible for extensive worldwide crop losses [3,39], other *Fusarium* species have demonstrated high incidence, causing severe and mild disease in tomato. *F. solani*, *F. equiseti*, *F. proliferatum*, and *F. verticillioides*, together with *F. oxysporum*, were already associated with rot diseases of tomato [12,13,16]. Akbar et al. [14,15] reported the prevalence of *F. equiseti* in tomato plants, which has been increasingly associated with many wilt diseases, followed by *F. graminearum*, *F. solani*, and *F. acuminatum*. These findings reinforce the idea of the importance of a single assay to target all the *Fusarium* species that might affect tomato.

Unexpectedly, our results allowed the detection and quantification of *Fusarium* spp. in several tomato plants from a commercial nursery grown under controlled conditions. Although these plant varieties present resistance to FOL race 3, these results are most likely explained by the lack of sanitation in the nurseries with the use of infected seeds and/or substrates, as early pathogen detection is crucial to assure the health status of the commercialized plant material [50]. The established assay might be used for the preventive detection of *Fusarium* diseases in plant material from nurseries, prior to transplantation into production fields.

5. Conclusions

Overall, the described *Fusarium* spp. qPCR TaqMan[®] assay targets a broad range of *Fusarium* species that cause extensive worldwide crop losses of tomato plants, allowing the discrimination of *Fusarium* spp. gDNA from other tomato pathogenic fungi. The assay revealed a high specificity, a high reliability, and a reproducible estimation of *Fusarium* accumulation in infected plants. The high sensitivity makes this methodology an efficient tool for the early diagnosis of the diseases, limiting the use of chemical defence strategies while providing an additional tool in the screening for resistant plants.

Supplementary Materials: The following supporting information can be downloaded at: <https://www.mdpi.com/article/10.3390/biology12020268/s1>, Figure S1: List of sequences of ITS region of *Fusarium* spp. and other fungal species affecting tomato in the Mediterranean Basin used on the phylogenetic analysis; Figure S2: Sensitivity and linearity of the *Fusarium* spp.-specific qPCR assay. Standard curve constructed based on Cq values obtained from the log of two-fold dilution series of gDNA from *F. oxysporum* f. sp. *radicis-lycopersici*. gDNA concentrations vary from 100 ng to 3.05×10^{-3} ng; Table S1: Quantification of *Fusarium* spp. in tomato field and control samples using the *Fusarium* spp.-specific qPCR assay. Data are expressed as quantification cycle (Cq) values.

Author Contributions: Conceptualization, M.D.C. and M.d.R.F.; methodology, M.D.C. and M.P.; formal analysis, M.D.C. and C.C.; investigation, M.D.C., C.V. and M.d.R.F.; resources, M.D.C., M.P., J.A.R., P.M. and A.A.; writing—original draft preparation, M.D.C.; writing—review and editing, C.V., M.P., J.A.R., C.C., P.M., A.A. and M.d.R.F.; funding acquisition, C.V. and M.d.R.F. All authors have read and agreed to the published version of the manuscript.

Funding: This work was funded by the projects “Development of a new virus-based vector to control TSWV in tomato plants” with the references ALT20-03-0145-FEDER-028266 and PTDC/ASPPLA/28266/2017, and “Control of olive anthracnose through gene silencing and gene expression using a plant virus vector” with the references ALT20-03-0145-FEDER-028263 and PTDC/ASPPLA/28263/2017, co-financed by the European Union through the European Regional Development Fund, under the ALENTEJO 2020 (Regional Operational Program of the Alentejo), ALGARVE 2020 (Regional Operational Program of the Algarve) and through the Foundation for Science and Technology (FCT), in its national component. Joana Ribeiro and Mariana Patanita are supported by Portuguese National Funds through FCT/MCTES under the PhD scholarships (2022.13638.BD and SFRH/BD/145321/2019, respectively). This work is also funded by National Funds through FCT under the Project UIDB/05183/2020.

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

Data Availability Statement: Data sharing not applicable.

Acknowledgments: The authors wish to thank to Sociedade Agrícola Pinhal das Teixeiras for providing the plant material used in the experiments.

Conflicts of Interest: The authors declare no conflict of interest.

References

1. Panno, S.; Davino, S.; Caruso, A.G.; Bertacca, S.; Crnogorac, A.; Mandi, A. A Review of the Most Common and Economically Important Diseases That Undermine the Cultivation of Tomato Crop in the Mediterranean Basin. *Agronomy* **2021**, *11*, 2188. [CrossRef]
2. Jones, J.B.; Zitter, T.A.; Momol, T.M.; Miller, S.A. *Compendium of Tomato Diseases and Pests, Second Edition*; The American Phytopathological Society: St. Paul, MN, USA, 2016; ISBN 978-0-89054-434-1.
3. Srinivas, C.; Nirmala Devi, D.; Narasimha Murthy, K.; Mohan, C.D.; Lakshmeesha, T.R.; Singh, B.P.; Kalagatur, N.K.; Niranjana, S.R.; Hashem, A.; Alqarawi, A.A.; et al. *Fusarium oxysporum* f. sp. *lycopersici* causal agent of vascular wilt disease of tomato: Biology to diversity—A review. *Saudi J. Biol. Sci.* **2019**, *26*, 1315–1324. [CrossRef] [PubMed]
4. Campos, M.D.; Félix, M.R.; Patanita, M.; Materatski, P.; Varanda, C. High throughput sequencing unravels tomato-pathogen interactions towards a sustainable plant breeding. *Hortic. Res.* **2021**, *8*, 171. [CrossRef] [PubMed]
5. Bodah, E.T. Root Rot Diseases in Plants: A Review of Common Causal Agents and Management Strategies. *Agric. Res. Technol. Open Access J.* **2017**, *5*, 555661. [CrossRef]
6. Yadeta, K.A.; Thomma, B.P.H.J. The xylem as battleground for plant hosts and vascular wilt pathogens. *Front. Plant Sci.* **2013**, *4*, 97. [CrossRef] [PubMed]

7. De la Lastra, E.; Basallote-Ureba, M.J.; De los Santos, B.; Miranda, L.; Vela-Delgado, M.D.; Capote, N. A TaqMan real-time polymerase chain reaction assay for accurate detection and quantification of *Fusarium solani* in strawberry plants and soil. *Sci. Hortic.* **2018**, *237*, 128–134. [\[CrossRef\]](#)
8. Ji, X.; Deng, T.; Xiao, Y.; Jin, C.; Lyu, W.; Wu, Z.; Wang, W.; Wang, X.; He, Q.; Yang, H. Emerging *Alternaria* and *Fusarium* mycotoxins in tomatoes and derived tomato products from the China market: Occurrence, methods of determination, and risk evaluation. *Food Control* **2023**, *145*, 109464. [\[CrossRef\]](#)
9. Maschietto, V.; Colombi, C.; Pirona, R.; Pea, G.; Strozzi, F.; Marocco, A.; Rossini, L.; Lanubile, A. QTL mapping and candidate genes for resistance to *Fusarium* ear rot and fumonisin contamination in maize. *BMC Plant Biol.* **2017**, *17*, 20. [\[CrossRef\]](#)
10. Debbi, A.; Bouregghda, H.; Monte, E.; Hermosa, R. Distribution and genetic variability of *Fusarium oxysporum* associated with tomato diseases in Algeria and a biocontrol strategy with indigenous *Trichoderma* spp. *Front. Microbiol.* **2018**, *9*, 282. [\[CrossRef\]](#)
11. McGovern, R.J. Management of tomato diseases caused by *Fusarium oxysporum*. *Crop Prot.* **2015**, *73*, 78–92. [\[CrossRef\]](#)
12. Rozlianal, F.S.; Sariah, M. Characterization of Malaysian Isolates of *Fusarium* from Tomato and Pathogenicity Testing. *Res. J. Microbiol.* **2006**, *1*, 266–272.
13. Murad, N.B.A.; Kusai, N.A.; Zainudin, N.A.I.M. Identification and diversity of *Fusarium* species isolated from tomato fruits. *J. Plant Prot. Res.* **2016**, *56*, 231–236. [\[CrossRef\]](#)
14. Akbar, A.; Hussain, S.; Ali, G.S. Germplasm Evaluation of Tomato for Resistance to the Emerging Wilt Pathogen *Fusarium equiseti*. *J. Agric. Stud.* **2018**, *5*, 174. [\[CrossRef\]](#)
15. Akbar, A.; Hussain, S.; Ullah, K.; Fahim, M.; Ali, G.S. Detection, virulence and genetic diversity of *Fusarium* species infecting tomato in Northern Pakistan. *PLoS ONE* **2018**, *13*, e0203613. [\[CrossRef\]](#)
16. Ribeiro, J.A.; Albuquerque, A.; Materatski, P.; Patanita, M.; Varanda, C.M.R.; Félix, M.R.; Campos, M.D. Tomato Response to *Fusarium* spp. Infection under Field Conditions: Study of Potential Genes Involved. *Horticulturae* **2022**, *8*, 433. [\[CrossRef\]](#)
17. Patanita, M.; Campos, M.D.; Félix, M.D.R.; Carvalho, M.; Brito, I. Effect of tillage system and cover crop on maize mycorrhization and presence of *Magnaportheopsis maydis*. *Biology* **2020**, *9*, 46. [\[CrossRef\]](#)
18. Le, K.D.; Kim, J.; Yu, N.H.; Kim, B.; Lee, C.W.; Kim, J.C. Biological Control of Tomato Bacterial Wilt, Kimchi Cabbage Soft Rot, and Red Pepper Bacterial Leaf Spot Using *Paenibacillus elgii* JCK-5075. *Front. Plant Sci.* **2020**, *11*, 775. [\[CrossRef\]](#)
19. Malik, M.S.; Haider, S.; Rehman, A.; Rehman, S.U.; Jamil, M.; Naz, I.; Anees, M. Biological control of fungal pathogens of tomato (*Lycopersicon esculentum*) by chitinolytic bacterial strains. *J. Basic Microbiol.* **2022**, *62*, 48–62. [\[CrossRef\]](#)
20. de Almeida, A.B.; Concas, J.; Campos, M.D.; Materatski, P.; Varanda, C.; Patanita, M.; Murolo, S.; Romanazzi, G.; Félix, M.R. Endophytic fungi as potential biological control agents against grapevine trunk diseases in alentejo region. *Biology* **2020**, *9*, 420. [\[CrossRef\]](#)
21. Heo, A.Y.; Koo, Y.M.; Choi, H.W. Biological Control Activity of Plant Growth Promoting *Rhizobacteria Burkholderia contaminans* AY001 against Tomato *Fusarium* Wilt and Bacterial Speck Diseases. *Biology* **2022**, *11*, 619. [\[CrossRef\]](#)
22. Campos, M.D.; Zellama, M.S.; Varanda, C.; Materatski, P.; Peixe, A.; Chaouachi, M.; Félix, M.R. Establishment of a sensitive qPCR methodology for detection of the olive-infecting viruses in portuguese and tunisian orchards. *Front. Plant Sci.* **2019**, *10*, 694. [\[CrossRef\]](#) [\[PubMed\]](#)
23. Azevedo-Nogueira, F.; Gomes, S.; Lino, A.; Carvalho, T.; Martins-Lopes, P. Real-time PCR assay for *Colletotrichum acutatum* sensu stricto quantification in olive fruit samples. *Food Chem.* **2021**, *339*, 127858. [\[CrossRef\]](#) [\[PubMed\]](#)
24. Lukianova, A.A.; Evseev, P.V.; Stakheev, A.A.; Kotova, I.B.; Zavriev, S.K.; Ignatov, A.N.; Miroshnikov, K.A. Development of qpcr detection assay for potato pathogen *Pectobacterium atrosepticum* based on a unique target sequence. *Plants* **2021**, *10*, 355. [\[CrossRef\]](#) [\[PubMed\]](#)
25. Thomas, W.J.; Borland, T.G.; Bergl, D.D.; Claassen, B.J.; Flodquist, T.A.; Montgomery, A.S.; Rivedal, H.M.; Woodhall, J.; Ocamb, C.M.; Gent, D.H. A Quantitative PCR Assay for Detection and Quantification of *Fusarium sambucinum*. *Plant Dis.* **2022**, *106*, 2601–2606. [\[CrossRef\]](#)
26. Campos, M.D.; Patanita, M.; Campos, C.; Materatski, P.; Varanda, C.M.R.; Brito, I.; Félix, M.R. Detection and quantification of *Fusarium* spp. (*F. oxysporum*, *F. verticillioides*, *F. graminearum*) and *Magnaportheopsis maydis* in maize using real-time PCR targeting the ITS region. *Agronomy* **2019**, *9*, 45. [\[CrossRef\]](#)
27. Varanda, C.M.R.; Materatski, P.; Landum, M.; Campos, M.D.; Félix, M.R. Fungal communities associated with peacock and cercospora leaf spots in olive. *Plants* **2019**, *8*, 169. [\[CrossRef\]](#)
28. Doyle, J.J.; Doyle, J.L. A rapid DNA isolation procedure for small quantities of fresh leaf tissue. *Phytochem. Bull.* **1987**, *19*, 11–15.
29. Varanda, C.M.R.; Oliveira, M.; Materatski, P.; Landum, M.; Clara, M.I.E.; Félix, M.R. Fungal endophytic communities associated to the phyllosphere of grapevine cultivars under different types of management. *Fungal Biol.* **2016**, *120*, 1525–1536. [\[CrossRef\]](#)
30. Edgar, R.C. MUSCLE: A multiple sequence alignment method with reduced time and space complexity. *BMC Bioinform.* **2004**, *5*, 113. [\[CrossRef\]](#)
31. White, T.J.; Bruns, S.; Lee, S.; Taylor, J. Amplification and direct sequencing of fungal ribosomal RNA genes for phylogenetics. *PCR Protoc. A Guide Methods Appl.* **1990**, *18*, 315–322.
32. Tamura, K.; Stecher, G.; Kumar, S. MEGA11: Molecular Evolutionary Genetics Analysis Version 11. *Mol. Biol. Evol.* **2021**, *38*, 3022–3027. [\[CrossRef\]](#)
33. Campos, M.D.; Valadas, V.; Campos, C.; Morello, L.; Braglia, L.; Breviario, D.; Cardoso, H.G. A TaqMan real-time PCR method based on alternative oxidase genes for detection of plant species in animal feed samples. *PLoS ONE* **2018**, *13*, e0190668. [\[CrossRef\]](#)

34. Materatski, P.; Varanda, C.; Carvalho, T.; Dias, A.B.; Campos, M.D.; Gomes, L.; Nobre, T.; Rei, F.; Félix, M.R. Effect of long-term fungicide applications on virulence and diversity of *Colletotrichum* spp. Associated to olive anthracnose. *Plants* **2019**, *8*, 311. [\[CrossRef\]](#)
35. Ruijter, J.M.; Ramakers, C.; Hoogaars, W.M.H.; Karlen, Y.; Bakker, O.; van den hoff, M.J.B.; Moorman, A.F.M. Amplification efficiency: Linking baseline and bias in the analysis of quantitative PCR data. *Nucleic Acids Res.* **2009**, *37*, e45. [\[CrossRef\]](#)
36. Bustin, S.A.; Benes, V.; Garson, J.A.; Hellemans, J.; Huggett, J.; Kubista, M.; Mueller, R.; Nolan, T.; Pfaffl, M.W.; Shipley, G.L.; et al. The MIQE guidelines: Minimum information for publication of quantitative real-time PCR experiments. *Clin. Chem.* **2009**, *55*, 611–622. [\[CrossRef\]](#)
37. Kralik, P.; Ricchi, M. A basic guide to real time PCR in microbial diagnostics: Definitions, parameters, and everything. *Front. Microbiol.* **2017**, *8*, 108. [\[CrossRef\]](#)
38. Hariharan, G.; Prasannath, K. Recent Advances in Molecular Diagnostics of Fungal Plant Pathogens: A Mini Review. *Front. Cell. Infect. Microbiol.* **2021**, *10*, 600234. [\[CrossRef\]](#)
39. Validov, S.Z.; Kamilova, F.D.; Lugtenberg, B.J.J. Monitoring of pathogenic and non-pathogenic *Fusarium oxysporum* strains during tomato plant infection. *Microb. Biotechnol.* **2011**, *4*, 82–88. [\[CrossRef\]](#)
40. Johnson, G.; Nolan, T.; Bustin, S.A. Real-time quantitative PCR, pathogen detection and MIQE. *Methods Mol. Biol.* **2013**, *943*, 1–16. [\[CrossRef\]](#)
41. Coats, K.; Debauw, A.; Lakshman, D.K.; Roberts, D.P.; Ismaiel, A.; Chastagner, G. Detection and Molecular Phylogenetic-Morphometric Characterization of *Rhizoctonia tuliparum*, Causal Agent of Gray Bulb Rot of Tulips and Bulbous Iris. *J. Fungi* **2022**, *8*, 163. [\[CrossRef\]](#)
42. Pavón, M.A.; González, I.; Martín, R.; García Lacarra, T. ITS-based detection and quantification of *Alternaria* spp. in raw and processed vegetables by real-time quantitative PCR. *Food Microbiol.* **2012**, *32*, 165–171. [\[CrossRef\]](#)
43. Lofgren, L.A.; Uehling, J.K.; Branco, S.; Bruns, T.D.; Martin, F.; Kennedy, P.G. Genome-based estimates of fungal rDNA copy number variation across phylogenetic scales and ecological lifestyles. *Mol. Ecol.* **2019**, *28*, 721–730. [\[CrossRef\]](#)
44. Lavrinienko, A.; Jernfors, T.; Koskimäki, J.J.; Pirttilä, A.M.; Watts, P.C. Does Intraspecific Variation in rDNA Copy Number Affect Analysis of Microbial Communities? *Trends Microbiol.* **2021**, *29*, 19–27. [\[CrossRef\]](#)
45. Kulik, T.; Jestoi, M.; Okorski, A. Development of TaqMan assays for the quantitative detection of *Fusarium avenaceum*/*Fusarium tricinctum* and *Fusarium poae* esyn1 genotypes from cereal grain. *FEMS Microbiol. Lett.* **2011**, *314*, 49–56. [\[CrossRef\]](#)
46. Bhagat, N.; Magotra, S.; Gupta, R.; Sharma, S.; Verma, S. Invasion and Colonization of Pathogenic *Fusarium oxysporum* R1 in *Crocus sativus* L. during Corm Rot Disease Progression. *J. Fungi* **2022**, *8*, 1246. [\[CrossRef\]](#)
47. Broeders, S.; Huber, I.; Grohmann, L.; Berben, G.; Taverniers, I.; Mazzara, M.; Roosens, N.; Morisset, D. Guidelines for validation of qualitative real-time PCR methods. *Trends Food Sci. Technol.* **2014**, *37*, 115–126. [\[CrossRef\]](#)
48. Perea-Domínguez, X.P.; Hernández-Gastelum, L.Z.; Olivas-Olguin, H.R.; Espinosa-Alonso, L.G.; Valdez-Morales, M.; Medina-Godoy, S. Phenolic composition of tomato varieties and an industrial tomato by-product: Free, conjugated and bound phenolics and antioxidant activity. *J. Food Sci. Technol.* **2018**, *55*, 3453–3461. [\[CrossRef\]](#)
49. Schrader, C.; Schielke, A.; Ellerbroek, L.; John, R. PCR inhibitors—occurrence, properties and removal. *J. Appl. Microbiol.* **2012**, *113*, 1014–1026. [\[CrossRef\]](#)
50. Azevedo-Nogueira, F.; Rego, C.; Gonçalves, H.M.R.; Fortes, A.M.; Gramaje, D.; Martins-Lopes, P. The road to molecular identification and detection of fungal grapevine trunk diseases. *Front. Plant Sci.* **2022**, *13*, 960289. [\[CrossRef\]](#)

Disclaimer/Publisher’s Note: The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.